



CHS-PU34


Phototherapy Unit




Cool comfortable lighting assures less fatigued eyes. Adjustable mobile stand designed for balanced stability & transporting. More effective phototherapy by 4pcs of fluorescent lamp.




• With Incubator



• Tilting



• Up-down



FL. 20BW-NU Lamp x 4ea
Total Energy : 10,440cal/h
400-500nm Energy : 5,600 cal/h

Specifications

Dimensions	620mm(W) x 540mm(D) 500mm(W) x 565mm(D) 1,150mm ~ 1,700mm(H)
Lamp	Non Ultra-violet ray Lamp 4pcs
Power source	AC 220V/110V, 50/60Hz
Power consumption	90W
Hour meter	0 ~ 999Hr
Caster	2"Swivel Caster 4pcs
Housing tilt	± 16°
Net weight	25kg

Standard Accessory Eye mask : 1pc / Curtain : 1pc