JW Medical

PRODUCT LIST

Jump to the World

Brighten The World
With JW Medical Technology

JW Values –"Respect for Life" and "the Pioneering Spirit"
JW Medical aims to lead advancement in healthcare services based on its founding principles of respect for life, and the pioneering spirit. Our purpose is to advance medical welfare through our own developments and cooperation with other leading global medical companies.

Operating Lights



Infant Incubator



OB/GYN Table



Operating Table



Digital X-ray



JW Medical

www.jw-medical.co.kr

Product List

02 Operating Room Equipment

04 Infant Care System

05 OB/GYN Equipment

06 Diagnostic Imaging Equipment

07 Sterilizer Equipment

08 Other Products

OPERATING ROOM EQUIPMENT



Honeylux LED RK

- Max illumination 160,000/ 120,000 Lux
- 10 Steps of dimming range
- 4,300K of color temperature
- >94Ra of general color rendering



Dialux Series

- Max illumination 150,000/ 80,000 Lux
- Hallogen type operating light
- 4,300K of color temperature
- 10 Steps of light adjusting mode



Honeylux LED

- Max illumination 160,000/ 120,000 Lux
- 10 Steps of dimming range
- 4,300K of color temperature
- >94Ra of general color rendering



Honeylux LED ECO

- Max illumination 130,000/ 90,000 Lux
- Dimming by IR remote controller (Option)
- >90Ra of general color rendering
- 5 Steps of dimming range
- 4,000K of color temperature



JW-T7000

- Wired & Wireless remote controller (Separate backup switch)
- Loading capacity 270kg in all positions
- Longitudinal shift 400mm
- Minimum height of 600mm (Base height 150mm)
- Auto locking system with 4 swivel casters



JW-T2000

- Remote control with simple user interface
- Loading capacity of 180Kg in all positions
- Manual longitudinal shift 300mm
- Minimum height of 605mm
- 50mm memory foam mattress



CHS-2500

- Longitudinal shift 300mm
- Height adjustment from 680 to 940mm
- C-ARM photography available to use
- X-Ray film cassette plate
- Wired remote controller



CHS-1500

- Electro-hydraulic multi-purpose operation table
- Height adjustment from 690 to 1,000mm
- C-ARM photography available to use
- X-Ray film cassette plate



CHS-790II

- Manual type operating table
- Height adjustment from 660 to 920mm
- Table-top rotating 360°



CHS-790

- Manual type operating table
- Height adjustment from 750 to 1000mm
- Table-top rotating 360°

INFANT CARE SYSTEM



JW-i3000

- 10.1" Full touch screen control
- Dual skin sensors
- Casette type of water reservoir
- Rail system for accessories
- SpO2 function (Masimo)



CHS-i1000

- Double wall type
- Air wall(curtain) function
- Air/skin temperature sensor
- Humidification function
- X-ray cassette tray with weight scale



CBW-1100

- Phototherapy lamp
- X-ray cassette tray
- Tilting mattress
- Head part 90° to either side for X-ray
- Alarm function for saftey



JW-PU1000

- 2 White LED/ 8 blue LED
- Height adjustment from 1200 to 1700(mm)
- Housing tilting ±35°
- Spectral Irradiance 60uW/cm2/nm at 40cm
- Set time from 30 to 240 min



CHS-PU34

- Non ultra-violet ray lamp 4pcs
- Housing tilting 16° up/down 1,150~1,700mm
- 400-500nm Energy: 5,600cal/h

OB/GYN EQUIPMENT



JW-G2000

- Hand remote controller
- Wired&wirelss foot switch system
- Automatic footrest
- Emergency battery(option)



CHS-E1000

- Gynecological examining table
- Wireless foot switch system
- Adjustable height 480~1,015mm
- Auto-position function



CHS-EX820

- Gynecological examining table
- Positioning by foot-switch
- Adjustable height 460~1,015mm
- Auto-position function



CHS-SMT204D


- Obstetric delivery table
- Linear actuator motor cylinder
- Height adjustment from 590 to 890mm
- Positioning by hand switch
- Memory function for positioning




CHS-E80

- Manual type obstetric delivery table
- 360° Horizontally rotating table top
- Height adjustment from 700 to 820mm
- Trendelenburg from -25° to 5°

DIAGNOSTIC IMAGING EQUIPMENT





VIDIX II

- Direct type dual detectors
- Auto tracking, synchronized, auto collimation function
- 7" LCD touch screen monitor
- Smooth 6-way motion table(2300mm*960mm)
- Tilting bucky stand



VIDIX B

- Direct type dual detectors
- 7" LCD touch screen monitor
- Simple design of bucky stand and table



VIDIX S

- Direct type single detector
- Detector smooth swing motion tilting 90°
- Mobile table
- Remote fuction control & foot switch

STERILIZER EQUIPMENT



CHS-ST045

- Steam sterilizaer (45 Liters)
- Squre type of chamber
- Self-diagnostic controlling total 13 errors conditions displayed
- LCD display
- Print function for option



CHS-ST065

- Steam sterilizaer (65 Liters)
- Squre type of chamber
- Self-diagnostic controlling total 13 errors conditions displayed
- LCD display
- Print function for option



CHW-EO55

- E.O gas sterilizaer (55 Liters)
- Round type of chamber
- Control & monitor display



CHS-AC80

- Steam sterilizaer (82 Liters)
- Safety door made from stainless steel plate
- Squre type of chamber
- Fully automatic control



CHS-AC250

- Steam sterilizaer (250 Liters)
- Safety door made from stainless steel plate
- Squre type of chamber
- Print function for option
- Fully automatic control

OTHER PRODUCTS



CHS-MATIC Series

- Multipurpose automatic scrub station
- Sterilized water (99.6%) supply - 2 U.V. water sterilizing lamps
- Hot water supply for four seasons
- Rust-proof materials-all stainless steel (SUS304)
- Dual filtering system



CHS-708

- Appliance for surgical operations and examinations
- 2 Suction bottles(made of P.C)
- Piston type oil-less pump
- Easy mobility by 4 castors
- 740mmHg ± 10mmHg of max. vacuum pressure



CHS-EV

- Appliance for obstetric and gynecological operations.
- 2 Suction bottles(made of P.C)
- Piston type oil-less pump
- Easy mobility by 4 castors
- 740mmHg ± 10mmHg of max. vacuum pressure

JW Medical

Jump to the World



Doctors and patients will be satisfied when using the high quality products from JW Medical

High-tech Medical Devices by JW Medical

JW Medical has developed high technology medical equipment through continuous study in the medical field. JW Medical helps to build a bright healthy society through constant development of their Infant Care, Operating Room and Digital X-ray systems.

Operating Lights



Infant Incubator



OB/GYN Table



Operating Table



Digital X-ray



www.jw-medical.co.kr



JW Medical

JW MEDICAL CORPORATION

HEAD OFFICE
JW Tower, 2477 Nambusunhwan-ro,
Seocho-gu, Seoul, Korea
Tel +82-2-2109-7846 Fax +82-2-852-0565
E-mail export@jw-medical.co.kr
www.jw-medical.co.kr

FACTORY
73, Chungjusandan 1-ro, Chungju-si,
Chungcheongbuk-do,
Republic of Korea (380-250)
Tel +82-43-855-7441